

### Shanghai Pigeon Office

📍 Add.: Room 2008, NO.8 Xianxia Road,  
Changning, Shanghai, China  
☎ +86 21 68880119  
☎ +86 21 68880117  
✉ pigeon@pigeondental.com

### Shanghai Pigeon Factory

📍 Add.: NO.233 Fengyong Road, Fengxian,  
Shanghai, China, 201414  
☎ +86 21 57561100  
☎ +86 21 57568302  
✉ pigeon@pigeondental.com

### Pigeon Japan Branch

📍 (203)1-35-14, Kitanaka, Hirakata City,  
Osaka, Japan  
☎ +81 72 835 7479  
☎ +81 72 835 7489  
✉ spj@pigeondental.com

pigeon@pigeondental.com

[www.pigeondental.com](http://www.pigeondental.com)



# SP-3 SUMMIT

## 3-Layer Composite Teeth

Product Catalog & Mould Chart

2016

[www.pigeondental.com](http://www.pigeondental.com)

# SP-3 SUMMIT Composite Teeth

The technology behind the smile:

Unique UNDER CUT technology (PATENT)

## Unique Under Cut Technology

- + The UNDER CUT design at the back of the teeth secures the teeth retention
- + Greatly improve the physical combining power of teeth and denture base
- + Makes a further integrate physical combination of the teeth and denture base besides the chemical covalent binding
- + Solve the difficulties of instable integration between the teeth and denture base in flexible denture making, achieve an obvious advantage in teeth fixing

## Excellent materials

- + Use Japan advanced composite materials and greatly improve the physical hardness of the teeth
- + With Outstanding Anti-dyeing property


## CAD / CAM Technology

- + Assures 100% accuracy of teeth dimension
- + Assures the accuracy of anatomical structure

## Precise Production Process

- + Without bubbles
- + Perfect smooth edges

Vickers Hardness ( Unit : HV0.3/15 )




# SP-3 SUMMIT

## 3-Layer Composite Teeth Mould Chart




# SP-3 SUMMIT

## 3-Layer Composite Teeth Mould Chart




# SP-3 SUMMIT POSTERIOR

## 3-Layer Composite Teeth Posterior Mould Chart



# SP-3 SUMMIT POSTERIOR

## 3-Layer Composite Teeth Posterior Mould Chart



### Shanghai Pigeon Office

📍 Add.: Room 2008, NO.8 Xianxia Road,  
Changning, Shanghai, China  
☎ +86 21 68880119  
📠 +86 21 68880117  
✉ pigeon@pigeondental.com

### Shanghai Pigeon Factory

📍 Add.: NO.233 Fengyong Road,Fengxian,  
Shanghai, China, 201414  
☎ +86 21 57561100  
📠 +86 21 57568302  
✉ pigeon@pigeondental.com

### Pigeon Japan Branch

📍 (203)1-35-14, Kitanaka, Hirakata City,  
Osaka, Japan  
☎ +81 72 835 7479  
📠 +81 72 835 7489  
✉ spj@pigeondental.com

pigeon@pigeondental.com

[www.pigeondental.com](http://www.pigeondental.com)



# SP-3 MARUJI 3-Layer Acrylic Teeth

Product Catalog & Mould Chart

2016

[www.pigeondental.com](http://www.pigeondental.com)

## SP-3 MALUJI 3-Layer Acrylic Teeth

The technology behind the smile:

Unique UNDER CUT technology (PATENT)

### Unique Under Cut Technology

- + The UNDER CUT design at the back of the teeth secures the teeth retention
- + Greatly improve the physical combining power of teeth and denture base
- + Makes a further integrate physical combination of the teeth and denture base besides the chemical covalent binding
- + Solve the difficulties of instable integration between the teeth and denture base in flexible denture making, achieve an obvious advantage in teeth fixing

### Excellent materials

- + Use Japan advanced synthetic resin materials and greatly improve the physical hardness of the teeth
- + With Outstanding Anti-dyeing property

### CAD / CAM Technology

- + Assures 100% accuracy of teeth dimension
- + Assures the accuracy of anatomical structure

### Precise Production Process

- + Without bubbles
- + Perfect smooth edges

# SP-3 MARUJI

## 3-Layer Acrylic Teeth Mould Chart




# SP-3 MARUJI

## 3-Layer Acrylic Teeth Mould Chart



# SP-3 MARUJI Posterior

## 3-Layer Acrylic Teeth Posterior Mould Chart



# SP-3 MARUJI Posterior

## 3-Layer Acrylic Teeth Posterior Mould Chart





### Shanghai Pigeon Office

📍 Add.: Room 2008, NO.8 Xianxia Road,  
Changning, Shanghai, China  
☎ +86 21 68880119  
☎ +86 21 68880117  
✉ pigeon@pigeondental.com

### Shanghai Pigeon Factory

📍 Add.: NO.233 Fengyong Road,Fengxian,  
Shanghai, China, 201414  
☎ +86 21 57561100  
☎ +86 21 57568302  
✉ pigeon@pigeondental.com

### Pigeon Japan Branch

📍 (203)1-35-14, Kitanaka, Hirakata City,  
Osaka, Japan  
☎ +81 72 835 7479  
☎ +81 72 835 7489  
✉ spj@pigeondental.com

pigeon@pigeondental.com

www.pigeondental.com



# SP-2 LUCK

## 2-Layer Acrylic Teeth

Product Catalog & Mould Chart

2016

www.pigeondental.com

# SP-2 LUCK

## 2-Layer Acrylic Teeth Mould Chart



# SP-2 LUCK

## 2-Layer Acrylic Teeth Mould Chart



# SP-2 LUCK POSTERIOR

## 2-Layer Acrylic Teeth Posterior Mould Chart



# SP-2 LUCK POSTERIOR

## 2-Layer Acrylic Teeth Posterior Mould Chart

